


# Monteringsanvisning H-Sele / H-Sele Mounting Description H-Belt bar / H-Belt H 型保护带及支架安装



1

## ENGLISH

Remove the push handles, and pull the H-Belt bar onto them (Picture 1.) Fix the push handles in required position. 拆去手推把,套入 H 型保护带支架(图 1)并调节至合适位置

The H-Belt bar is telescopic, so pull it out to required width. On both plastic reels there is a setscrew. Tighten these screws with a 2,5mm Allen key. Fix the H-Belt bar tightening the setscrews on each side using a 2.5mm Allen key (Picture 2). The height should be level with the shoulders of the user.

H 型保护带支架带有伸缩功能, 可将其调到所需的宽度  
在两个侧塑料卷轮上均带有紧定螺丝, 可使用 2.5MM 内六角将其可松开或扭紧(图 2).  
支架高度应与使用者肩膀平行.



1

2/3



Thread the belt through the reels and lock the belt by pulling it through the belt clamp (Picture 3). Adjust to the required length of belt. 将带子穿过支架上的塑料卷轮后, 再穿入(塑料)带夹中并拉紧(图 3),调节带子长度.

4

The lower part, the hip belt, is mounted as a standard hip belt. Pull the belt through the hole in the hip belt bracket (Picture 4).

低位处, 那段臀部带子安装方法与标准臀部保护带安装方法相同. 将带子套入臀部处的带扣孔中(图 4)

4

Thread the belt back through the enclosed belt clamp (Picture 5).


将带子穿过带夹

5

Fix the hip belt bracket in the bacmost hole of the backrest hinge using the enclosed screws and nuts. Use two 13mm open-end spanners.

用所配螺丝和螺母将臀部带扣安装在靠背支架最后一个孔上.

工具: 2 把 13MM 开口扳手



5



