



Product Service

CERTIFICATE

No. Z1 12 03 22430 034

Holder of Certificate: Alu Rehab ASBedriftsveien 23
4353 Klepp Stasjon
NORWAY**Certification Mark:****Product:****Standard wheelchairs
Manual wheelchair**

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.:

713001019

Date, 2012-03-19

(Torsten Zimmer)

Page 1 of 2

**CERTIFICATE**

No. Z1 12 03 22430 034

Model(s): **Netti 4U Comfort CE PLUS**

Parameters:

total length incl. footrests:	1130 mm
total length without footrests:	930 mm
width:	530 - 680 mm
height:	1100 mm
seat width:	350, 400, 450 and 500 mm
seat depth:	425 - 500 mm
seat height:	465 - 500 mm
backrest height:	500 mm
seat plane angle:	-5° - +20°
backrest angle:	92° - +137°
max. user weight:	135 kg
brake system:	knee lever brake (user brake)
back wheels:	
- size	24 x 13/8"
- tyres	accident proof
- adjustable	Yes
- fixation	Quick release
front wheels:	
- size	7 x 1 3/4", Ø 175 x 33 mm
- tyres	accident proof
- adjustable	Yes
- fixation	Quick release
frame material:	aluminium
armrest:	adjustable and removable
legrest:	height, angle adjustable and removable
anti tip device:	adjustable

Tested according to: DIN EN 12183:2009**Production Facility(ies):** 57370