



Product Service

CERTIFICATE

No. Z1 10 11 22430 024

Model(s): **Manual wheelchair type: Netti III Comfort**

Parameters:

total length	
incl. footrests:	1100 mm
total length	
without footrests:	900 mm
width:	590 mm – 740 mm
height:	1120 mm
seat width:	350, 380, 400, 430, 450 and 500 mm
seat depth:	400 mm - 500 mm
seat height:	400 mm - 480 mm
backrest height:	480 mm (600 mm with back extension)
seat plane angle:	- 9° - +16°
backrest angle:	86° - 133°
max. load:	145 kg
brakes:	knee lever brake (user brake), drum brake (assistant brake)
back wheels	
-size:	20 x 1", 22 x 1", 24 x 1" and 24 x 1 3/8"
-tyres:	solid rubber (accident proof)
-adjustable:	yes
-fixation:	Quick release
front wheels	
-size:	∅ 125 x 27, ∅ 150 x 27 and 175 x 33 mm
-tyres:	solid rubber (accident proof)
-adjustable:	yes
-fixation:	Quick release
frame material:	aluminium
armrest:	height adjustable and removable
legrest:	height adjustable, angle adjustable and removable
anti tip device:	adjustable

Tested according to: DIN EN 12183:2009

Production Facility(ies): 57370



Product Service

CERTIFICATE

No. Z1 10 11 22430 024

Holder of Certificate: Alu Rehab ASBedriftsveien 23
4353 Klepp Stasjon
NORWAY**Certification Mark:****Product:****Standard wheelchairs
Manual wheelchairs**

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.:

71375626

**Date,** 2010-11-24
(Torsten Zimmer)

Page 1 of 2